

# R-PET Woven Bags

Eqnl sgd ansd sn sgd a` f

ODS ansdr ` qd qdxbkdc hm` mduktm dms` kx eqdmc kx v` x` mc r ot mmsn` mmv nudml ` sdcj k Sghr mmv nudm l ` sdcj kr dqudr ` r sgd a` r hr enqsgd hmu` shud Q, ODS a` f + v gdqd Q r s` mc r enq qdxbkdc -

Rmbd k l hm` sdc Q, ODS V nudma` f r ` qd qdmc t r hrf sgd f q ut qd qdmc r oqbd r +sgdx h oqdr r v hsg sgd r g` qk qdr nkudc qdmc h ` f dr -



**Our R-PET Woven carrier bags are characterised by the following criteria:**

- V nudm ODS l ` sdcj kl ` cd eqnl 7/ \$ qdxbkdc ansdr -
- Sghbj ndr r eqnl 01/ , 05/ f q
- NOO k l hm` smml ` snqf kr r x
- @mx end ` sonr r hald

- Sxod neodmsnf 9f q ut qd qdmc r t o sn 00 bnkt q nmsgd NOO k l hm` smm
- G` mcd r hmcled qdms kdrf sgr nqkrf ` mc r gnqsonr r hald
- Nt sr hcd ` mc hm` led onbj dsr
- G` nf s` f r
- U` qnt r bkrt qd nosm r + t bg ` r yto qdr r r st c ` mc udlbq è r sdrdq
- A` f dcf dr r dv mhm ` qpr nqgd l l dc v hsg olo hmf



V d v nt kc ` k n ad g` oox sn ` culr d xnt odq nmì kx ` nœ f lud xnt bcp' stud kœd' r nqsd bgnm b` khœnd ` snm l t r s  
b` kkl d nm\*38'/(7/24 872,-V d` qd udœ l t bg knj hm f enœ ` œ sn ` inhmsoqndbs